

BODY AESTHETICS

Premium Liposuction Cannulas for Body Contouring



+923216126225



info@adasinstrumentsonline.com
adas.surgical@gmail.com



www.adasinstrumentsonline.com



Table Of Contents

- 1- Injector Liposuction Cannula Page: 2-3
- 2- Liposuction Accessories Page: 4
- 3- Cannula Hubs / Handle info Page: 5
- 4- General Suction / Harvesting Cannula Page 6-7
- 5- Infiltration Cannula Page: 8 to10
- 6- Accessories and handles Page: 11
- 7- Customize Sets Page:12 to17





Luer Lock



Super Luer Lock



Toomey

Standard Micro Injector



Two Hole Micro Injector



Blunt side-port tip glides effortlessly through soft tissue while protecting vascular structures — ideal for small volume fat injection.

Superficial Injector



Five graduated spiral ports create a fan-shaped microfat distribution for smooth, even intradermal injection results.

Spoon Tip Injector



Blunt spoon-shaped tip with a wide port for rapid, high-volume fat injection into larger treatment areas.

Toledo V Dissector



V-shaped blunt tip designed for precise superficial tissue dissection and smooth fat injection in a single instrument.

Little injector



Single small round port for ultra-precise microfat delivery in delicate areas like the eye region and lip border.

Sharp extractor / Injector



Blunt extractor / Injector



Dual-function tip designed for both precise adipose tissue extraction and smooth fat reinjection in a single instrument

Size	Micro Injector	Two Hole	Superficial	Spoon Tip	Toledo V Dissector	Little injector	Sharp extractor / Injector	Blunt extractor / Injector
0.7mm (22g) x 4cm	010704	020704	030704			060704		
0.7mm (22g) x 7cm	010707	020707	030707			060707		
0.9mm (20g) x 4cm	010904	020904	030904			060904		
0.9mm (20g) x 7cm	010907	020907	030907			060907		
1.0mm (19g) x 10cm	011010	021020	031030	041040		061060	071070	081080
1.0mm (19g) x 15cm	011015	021025	031035	041045		061065	071075	081085
1.2mm (18g) x 10cm	011210	021210	031210	041210		061210	071210	081210
1.2mm (18g) x 12cm	011212	021212	031212	041212		061212	071212	081212
1.2mm (18g) x 15cm	011215	021215	031215	041215		061215	071215	081215
1.4mm (17g) x 10cm	011410	021410	031410	041410		061410	071410	081410
1.4mm (17g) x 15cm	011415	021415	031415	041415		061415	071415	081415
1.4mm (17g) x 20cm	011420	021420	031420	041420		061420	071420	081420
1.6mm (16g) x 15cm	011615	021615		041615	051615	061615	071615	081615
1.6mm (16g) x 20cm	011620	021620		041620	051620	061620	071620	081620
1.6mm (16g) x 25cm	011625	021625		041625	051625	061625	071625	081625
2.0mm (12g) x 12cm				042042	052052		072072	082082
2.0mm (12g) x 15cm				042045	052055		072075	082085
2.0mm (12g) x 20cm				042020	052020		072020	082020
2.5mm (10g) x 15cm				042415	052415		072415	082415
2.5mm (10g) x 20cm				042420	052420		072420	082420
2.5mm (10g) x 25cm				042425	052425		072425	082425

Select Code of hub and Cannula
For order ie A-012015

NOTE: When You Need Curve Cannula

Just Add Curve Code With Product Code Like A-010704CB



Curved Injectors

These Curved Injectors are designed especially for fat transfer procedures to the body.

Here you can find most regular curved shape however you can draw your own shape and send us too we will made.

Follow are the Regular shape:

- 1: Angle 45 degree bent
- 2: Round curved injector with Concave Hole
- 3: Round Curved injector with Convex hole
- 3: S Shape curved Injector cannula





NANO Fat Filter

Art #	Size
A-507-01	Nano Filter
A-507-02	Filter 0.25
A-507-03	Filter 0.50
A-507-04	Filter 0.75

Fine-mesh filtration system designed to process and purify harvested fat into ultra-refined nano fat, removing fibrous tissue and impurities for smoother, more precise micro-injection results.

Lure to Lure (Anaerobic Transfers)



Airtight luer-to-luer connector enabling seamless, oxygen-free fat transfer between syringes, preserving fat cell viability for improved grafting survival rates.

Art #	Size	Art #	Size
A-500.01	1.0mm hole	A-500.06	2.4mm hole
A-500.02	1.2mm hole	A-500.07	2.5mm hole
A-500.03	1.4mm hole	A-500.08	3.0mm hole
A-500.04	1.5mm hole	A-500.09	0.8mm three hole
A-500.05	2.0mm hole	A-500.10	1mm three hole

Syringe Stopper Cap



Protective cap that seals the syringe tip securely, preventing contamination and maintaining sterility during fat storage and transfer procedure

Art #	Size
A-501.01	Aluminium
A-501.02	Steel

Suction Machine Adapter Attachment For Luer Lock and Super Luer Lock Cannula



Universal adapter attachment that connects standard luer lock and super luer lock cannulas directly to suction machines, ensuring a secure, airtight connection for consistent and uninterrupted suction performance during liposuction procedures.

Art #	Size
A-508.01	For luer lock cannula

I.V Spikes



Art #	Size
A-509.01	Steel
A-509.02	Aluminium fix

Sterile sharp-tip spike designed for quick and secure insertion into IV fluid bags, enabling smooth and uninterrupted fluid flow during surgical infiltration procedures.

Fact Injector Gun



Art #	Size
A-506-01	Sliver 10 / 20 CC
A-506-02	Black 10 / 20 CC
A-506-03	Sliver 50 / 60 CC
A-506-04	Black 50 / 60 CC

Precision-engineered fat injection gun for smooth, controlled fat transfer during grafting procedures, delivering consistent pressure and accurate volume output with every stroke

Syringe Stop Lock



Art #	Size
A-502-01	10 cc
A-502-02	20 cc
A-502-03	50 cc
A-502-04	60 cc

Secure locking mechanism that holds the syringe plunger in place, maintaining consistent pressure and preventing accidental fluid release during fat transfer procedures.

Female Luer conector (hose adapter)



Art #	Size
A-511.01	For luer lock cannula

Secure female luer connector that links the cannula or tubing to suction hoses, ensuring a leak-proof, airtight connection for smooth and uninterrupted fluid flow during procedures.

Holding device for Luer Lock syringes

Art #	Size
A-505.01	Stain
A-505.02	Gold Plated

Ergonomic holding device designed to firmly secure luer lock syringes in place during fat transfer procedures, providing stable one-handed control for precise and consistent injection delivery.



Automatic Taralock Fat Syringe



Art #	Size
A-513.01	10cc
A-513.02	20cc
A-513.03	60cc

Select Code of hub, Cannula, and Size For order ie A-012015

Luer Lock (Aluminum and Steel)



Luer Lock Hub A secure threaded connection system that firmly locks the cannula to the syringe, preventing accidental disconnection during procedures. The twist-and-lock mechanism ensures an airtight, leak-proof seal for consistent pressure control throughout fat harvesting and injection. Compatible with all standard syringes, offering surgeons a reliable and safe connection for uninterrupted workflow.

Super Luer Lock Hub An enhanced heavy-duty version of the standard luer lock, featuring a reinforced threading system for an even more secure and pressure-resistant connection. Ideal for high-volume procedures where maximum stability and zero fluid leakage are critical for precise, uninterrupted performance.

Super Luer Lock



Toomey 60cc Hub



Toomey 60cc Hub Large-capacity 60cc hub designed for high-volume fat harvesting and transfer procedures, reducing the need for frequent syringe changes. The wide-bore connection ensures smooth, unrestricted flow for faster and more efficient fat collection in large-scale body contouring cases.

Threaded Fitting Hub A precision-threaded connection designed for a firm and stable instrument attachment, ensuring consistent performance throughout the procedure. The secure thread design eliminates unwanted movement or rotation, providing surgeons full control and confidence during both extraction and injection.

Threaded Fitting



Machine One Piece Handle (Super Handle)



Super Handle A heavy-duty machine-compatible handle engineered for secure attachment to powered liposuction systems, delivering stable and vibration-free performance during high-speed procedures. The reinforced grip and robust connection ensure precise directional control, reducing surgeon fatigue in long and demanding surgical sessions. Available as a solid one-piece construction for maximum strength and durability.

Single-Port Classic Cannula

Code: 01



The most widely used cannula design in body contouring. Features one extraction port positioned opposite the thumb grip for maximum surgeon control and consistent fat removal.

Two Hole

Code: 02



Three Hole

Code: 03



Features two or three extraction ports aligned straight along the cannula shaft for smooth and balanced fat removal.

Mercedes

Code: 04



Double Mercedes

Code: 05



Triple Mercedes

Code: 06



Mercedes Cannula Series (Single, Double & Triple) Featuring circumferential ports in 3, 6, and 9-port configurations, this series delivers smooth multi-plane fat removal for every procedure size from precise contouring to high-volume extraction.

Pyramid Hole Cannula

Code: 07



Sharp pyramid-shaped tip with precision ports designed to penetrate and break up dense, fibrous fat deposits with ease.

Stevens Speed Style Cannula

Code: 08



Multiple small ports in a spiral configuration for advanced liposculpture, delivering fast and efficient fat removal with minimal tissue trauma.

Three Spiral Port

Code: 09



Six Spiral Port

Code: 10



Three Spiral Port HarvesterThree ports in a spiral layout for smooth, controlled fat harvesting across multiple tissue planes in a single pass.

Six Spiral Hole HarvesterSix spiral-arranged ports for high-volume fat collection, significantly reducing operative time in large-scale harvesting procedures.

Basket Cannula

Code: 11



Double Basket Cannula

Code: 12



Triple Basket Cannula

Code: 13



Featuring multi-port basket configurations in single, double, and triple arrangements, this series excels in areas with skin irregularities offering precise injection, suction, and tissue equalization across all procedure sizes.

Accelerator Cannula

Code: 14



Engineered for rapid, high-efficiency fat removal

Cobra Bibevel Style

Code: 15



For Dense, fibrous tissue with precision and control.

Cobra K

Code: 16



Effective for fat removal in dense fibrous tissue areas.

V Dissector

Code: 17



V-shaped cutting tip designed to release and break apart dense fibrous tissue in the superficial compartment with precision.

Becker Grater Round

Code: 18



Micro-perforated round tip designed for smooth, uniform fat removal with minimal tissue disruption.

Sturdy multi-tier stainless steel syringe stand featuring 17 organized holes to accommodate multiple syringe sizes from 10cc to 60cc. Designed to keep syringes securely upright and within easy reach during surgical procedures for a clean, efficient, and well-organized workspace.

Sryinge Stand

Art #	Size
A-509.10	Sryinge Stand



Gynomastia cannula

Code: 19



designed for targeted fat and glandular tissue removal in male chest contouring procedures.

Toledo V Dissector

Code: 20



Toledo V Dissector Double hole

Code: 21



for precise subcutaneous tissue dissection and smooth fat release in superficial treatment areas.

Cane Superficial Liposuction

Code: 22



controlled fat removal in superficial skin layers for smooth and refined body contouring results.

Triport Small Harvesting

Code: 23



Triport II

Code: 24



Three-port arrangement delivering enhanced fat collection capacity while maintaining precision control in confined treatment areas.

Helixed Tri-port

Code: 25



Three ports arranged in a helical pattern for smooth, evenly distributed fat extraction across multiple tissue planes with minimal resistance.

KOTZUR 9 Holes

Code: 26



Kotzur Cannula Nine strategically distributed ports including large oval tip openings for powerful initial suction, delivering effective and gentle fat removal in densely connected tissue areas such as gynecomastia treatment.

Spatula Single port

Code: 27



for precise subcutaneous tissue dissection and smooth fat release in superficial treatment areas.

One Central Two Lateral Ports

Code: 28



Two Lateral One Central Ports

Code: 29



High-volume fat removal with enhanced directional suction control.

Tonnard Harvester

Code: 30



Elevated sharp microports in a spiral arrangement for efficient fat cell harvesting on both forward and backward strokes.

Sorenson Harvester

Code: 31



Recessed sharp microports designed for gentle yet efficient fat harvesting in both push and pull directions.

Flared Mercedes

Code: 32



Multiple circumferential hole pattern along the shaft with a rounded tip for maximum fat intake and uniform extraction across all tissue planes.

Size Chart

Code	Size	Code	Size
2015	2.0mm (12g) x 15cm	4030	4.0mm (6g) x 30cm
2020	2.0mm (12g) x 20cm	4035	4.0mm (6g) x 35cm
2515	2.5mm (10g) x 15cm	5020	5.0mm (4g) x 20cm
2520	2.5mm (10g) x 20cm	5025	5.0mm (4g) x 25cm
2525	2.5mm (10g) x 25cm	5030	5.0mm (4g) x 30cm
3015	3.0mm (8g) x 15cm	5035	5.0mm (4g) x 35cm
3020	3.0mm (8g) x 20cm	6040	6.0mm (2g) x 40cm
3025	3.0mm (8g) x 25cm	6045	6.0mm (2g) x 45cm
4015	4.0mm (6g) x 15cm	8035	8.0mm (0g) x 35cm
4020	4.0mm (6g) x 20cm	8040	8.0mm (0g) x 40cm
4025	4.0mm (6g) x 25cm	8045	8.0mm (0g) x 45cm


Select Code of hub & Cannula & Size
For order ie A-012015

Our liposuction cannulas are manufactured to the highest standards, available in custom dimensions, sizes and configurations. We offer full private labeling and OEM options for clinics and distributors.

Infiltration Cannula For tumescent local anesthesia Procedures

All Tumescent Infiltrators featured are precision-designed for tumescent local anesthesia in liposuction procedures. Available in various hole patterns to ensure smooth, even and controlled fluid delivery throughout the subcutaneous tissue.

18 Hole **22** Hole **26** Hole



Infiltration Cannula

Code	Lenght	Ø	holes	Ø Hole
011520	15 cm/6"	2,0 mm	18	1,0 mm
011525	15 cm/6"	2,5 mm	18	1,2 mm
011530	15 cm/6"	3,0 mm	18	1,5 mm
012020	20 cm/8"	2,0 mm	22	1,0 mm
012025	20 cm/8"	2,5 mm	22	1,2 mm
012030	20 cm/8"	3,0 mm	22	1,5 mm
012040	20 cm/8"	4,0 mm	22	2,0 mm
012525	25 cm/10"	2,5 mm	26	1,2 mm
012530	25 cm/10"	3,0 mm	26	1,5 mm
012540	25 cm/10"	4,0 mm	26	2,0 mm
013030	30 cm/12"	3,0 mm	26	1,5 mm
013040	30 cm/12"	4,0 mm	26	2,0 mm
013050	30 cm/12"	5,0 mm	26	2,5 mm

Infiltration Cannula For tumescent local anesthesia Procedures


26 Hole



Tumescent Infiltrator Angled

Code	Lenght	Ø	holes	Ø Hole
032030	20 cm/8"	3,0 mm	26	1,5 mm
032040	20 cm/8"	4,0 mm	26	2,0 mm
032530	25 cm/10"	3,0 mm	26	1,5 mm
032540	25 cm/10"	4,0 mm	26	2,0 mm
033030	30 cm/12"	3,0 mm	26	1,5 mm
033040	30 cm/12"	4,0 mm	26	2,0 mm
033050	30 cm/12"	5,0 mm	26	2,5 mm

30 Hole



Code	Lenght	Ø	holes	Ø Hole
022030	20 cm/8"	3,0 mm	30	1,5 mm
022040	20 cm/8"	4,0 mm	30	2,0 mm
022530	25 cm/10"	3,0 mm	30	1,5 mm
022540	25 cm/10"	4,0 mm	30	2,0 mm
023030	30 cm/12"	3,0 mm	30	1,5 mm
023040	30 cm/12"	4,0 mm	30	2,0 mm
023050	30 cm/12"	5,0 mm	30	2,5 mm

8 Hole **10** Hole

For fine proportioning in facial and neck areas

Code	Lenght	Ø	holes	Ø Hole
041208	12 cm/4-1/2"	2,0 mm	8	1,0 mm
041508	15 cm/6"	2,0 mm	8	1,0 mm
041210	12 cm/4-1/2"	2,0 mm	10	1,0 mm
041510	15 cm/6"	2,0 mm	10	1,0 mm



Infiltration Cannula

Infiltration Cannula For tumescent local anesthesia Procedures



Bayonet-shaped design with offset angle provides superior access to confined spaces while keeping hands clear for optimal visibility and control

22 Hole

30 Hole

Code	Lenght	Ø	holes	Ø Hole
052030	20 cm/8"	3,0 mm	22	1,5 mm
052040	20 cm/8"	4,0 mm	22	2,0 mm
053030	30 cm/12"	3,0 mm	30	1,5 mm
053040	30 cm/12"	4,0 mm	30	2,0 mm
053050	30 cm/12"	5,0 mm	30	2,5 mm

Select Code of hub, and Cannula,
For order ie A-052030

Bayonet Infiltration Cannula

Custom Infiltration Cannulas - Built to Your Specifications

Share your requirements and we'll manufacture to your specifications. Fully customizable to meet your exact clinical needs - select from multiple gauges, lengths, tip designs, and side-port configurations for optimal fluid delivery, precise tissue access, and superior procedural control.



Decant Stand



Art #	Size
A-511.01	Luer to Luer
A-511.02	60cc to Luer
A-511.03	60cc to 60cc

Designed for the safe handling and movement of harvested fat tissue prior to reinjection, the Decant Stand utilizes gravity to naturally separate and decant the extracted material without manual intervention

Hose adapter



Art #	Size
A-512.01	10 cc
A-512.02	20 cc
A-512.03	60 cc

The Hose Adapter connects a syringe and cannula to an aspirating machine, allowing for efficient and controlled tissue aspiration. Compatible with both 20cc and 60cc syringes for versatile clinical use.

Automatic Taralock Fat Syringe



Art #	Size
A-513.01	10cc
A-513.02	20cc
A-513.03	60cc

The Automatic Taralock Fat Syringe is a precision instrument developed for manual liposuction and fat transfer procedures. It automatically generates negative pressure during aspiration, reducing surgeon effort while maintaining consistent suction. The built-in power injector enables smooth re-injection of harvested soft tissue through both large and small diameter cannulas with minimal force and maximum control.

Skin Punch Sharp Tip

Art #	Size
A-514.01	12cm
A-514.02	15cm
A-514.03	18cm



Handle with connection thread

The Threaded Cannula Handle is designed for use in liposuction procedures, featuring an internal threaded connection that securely fits into threaded liposuction cannulas. The handle provides a firm and stable grip while maintaining a sealed connection to the suction machine for consistent and controlled aspiration.



Art #	Size
A-504.01	Ø 13 mm
A-504.02	Ø 13 mm with hole

Handle with Lure connection thread

A compact dual-connection handle designed for liposuction procedures. One end features an internal thread that securely fits into threaded liposuction cannulas, while the opposite end accepts a luer lock syringe for precise, syringe-based aspiration and tissue transfer.



Art #	Size
A-503-01	with lure lock

FAT Collection Jar

A high-capacity sterile collection jar designed for efficient fat harvesting during liposuction procedures. The sealed canister maintains a closed system to preserve fat cell viability, making it suitable for both fat removal and fat grafting applications. Constructed from medical-grade glass with a durable lid, it is fully reusable and autoclavable for repeated clinical use.



Art #	Size
A-515.01	3000 ml

Facial Set AD-5000

A complete all-in-one set designed for fat harvesting and transfer procedures, combining precision instrumentation for both aspiration and re-injection. The set features ultra-fine super luer lock cannulas across multiple diameters and lengths to accommodate varying tissue depths, paired with anaerobic transfer connectors to maintain a closed sterile system. A titanium-coated Sorensen harvester ensures gentle and efficient fat collection, while the included nano fat set and fine facial cannula extend the kit's versatility to delicate contouring and micro-fat grafting applications.



Description	Qty
Canula 0.9 mm x 4cm - Super Luer Lock	1
Canula 1.2 mm x 10cm Super Luer Lock	1
Canula 1.4 mm x 15cm Super Luer Lock	1
Lure to Lure (Anaerobic Transfers) 1mm	1
Lure to Lure (Anaerobic Transfers) 2mm	1
Fine Facial and neck cannula 2.5mm x 20cm 10 hole	1
Sorensen Harvester 1.0 mm Holes - 20 cm - Titanium coated - Luer lock	1
Nano Fat set	1
Decant Stand Anerobic Transfer Toomey 60cc to Luer	1

Facial & Neck Liposuction Cannula Set AD-5010

A fully comprehensive set engineered for precise facial and neck liposuction and fat grafting procedures. The kit brings together a wide selection of specialized luer lock cannulas — including Mercedes tip, Spatula tip, Kotzur gynecomastia, tumescent infiltrator, V-Dissector, and multi-hole harvester designs — covering every stage from infiltration and aspiration to tissue harvesting and re-injection. Complemented by luer lock transfer adaptors, a dedicated handle, syringe stoppers, and a 60cc syringe stop lock, this set provides the surgeon with a complete, ready-to-use solution for refined and controlled body contouring of delicate facial and neck areas.

Description	Qty
3mm x 25cm KOTZUR gynecomasty 9 holes Cannula, lure lock	1
3mm x 25cm Mercedes tip lure lock liposuction cannula	1
3mm x 25cm Tumescent Infiltrator curved lure lock liposuction cannula	1
4mm x 25cm Three standard hole lure lock liposuction cannula	1
3mm x 25cm Spatula tip lure lock liposuction cannula	1
2mm X 25cm Harvester cannula with shevels hole lure lock	1
2.5mm x 20cm Harvester cannula with spikes hole Super lure l	1
2.5mm x 20cm two lateral hole lure lock liposuction cannula	1
4mm x 15cm V-Dissector Liposuction cannula lure lock	1
3mm x 15cm Harvester 6 small hole lure lock liposuction Cannula	1
2mm x 15cm reinjecting cannula lure lock	1
1mm x 5cm reinjecting cannula Super lure lock	1
1mm x 10cm reinjecting cannula lure lock	1
Lure lock to Lure lock transfer adaptor	2
Handle for luer lock	1
syringe Stopper	2
syringes Stop lock 60cc	1



All sets can be customized to suit your individual surgical needs. Cannulas and instruments may be added, removed, or replaced according to your preference. Please contact us for your personalized configuration.

Multipurpose Cannula Set Facial, Breast, Thigh & Buttock AD-5011

Designed to meet the demands of modern aesthetic surgery, this multipurpose set equips the surgeon to perform liposuction, fat harvesting, and precision fat grafting across multiple body areas in a single, cohesive kit. From delicate facial contouring to larger volume procedures of the breast, thigh, and buttock, every instrument has been carefully selected to ensure smooth workflow, consistent results, and maximum versatility in the operating theatre.



Description	Qty	Description	Qty
Lure lock to Lure lock transfer adaptor	2	2.5mm x 20cm two lateral hole lure lock liposuction cannula	1
Lure lock to Lure lock transfer adaptor	1	3.0mm x 20cm Blugerman I luer lock cannula	1
iv spike adapter	1	2.5mm X 20cm Harvester cannula with shevels hole lure lock	1
Handle for luer lock	1	2.1mm x 20cm Tumescant Infiltrator luer lock cannula	1
Nano Fat with 3 filter	1	3mm x 20cm Mercedes tip Super lure lock liposuction cannula	1
syringe Stopper	2	4mm x 20cm Spatula tip lure lock liposuction cannula	1
syringes Stop lock 10cc / 20cc / 50cm	3	3mm x 25cm Tumescant Infiltrator luer lock cannula	1
0.9mm x 5cm reinjecting cannula lure lock	1	4mm x 25cm Three standard hole lure lock liposuction cannula	1
0.7mm x 4cm reinjecting cannula lure lock	1	joseph skin hook double prong	1
2.0mm x 5cm reinjecting cannula lure lock	1	Iris Super cut scissors 11.5cm	1
1.2mm x 5cm reinjecting cannula lure lock	1	Neivert needle holder Tungsten Carbide 12.5cm	1
2.0mm x 15cm reinjecting cannula Super luer lock	1	Gorney Platysma scissors Tungsten Carbide 17cm	1
3mm x 15cm V?Dissector Liposuction cannula lure lock	1		

Multipurpose Cannula Set Facial, Breast, Thigh & Buttock AD-5012

A carefully assembled set built for surgeons who require versatility across multiple aesthetic procedures. Covering liposuction, fat harvesting, and precision fat grafting for the face, breast, thigh, and buttock, the set includes cannulas ranging from ultra-fine 0.7mm re-injection tips for delicate micro and nano fat grafting all the way to 4mm liposuction cannulas for larger volume aspiration. Diverse tip configurations including Mercedes, Spatula, Blugerman I, V-Dissector, Tumescant Infiltrator, and Harvester designs ensure the right instrument is always at hand for every anatomical zone and tissue depth. A fully integrated accessory kit including nano fat filtration, luer lock transfer adaptors, IV spike adapter, syringe stoppers, and stop locks makes this a complete theatre-ready solution.



Description	Qty
0.9mm x 5cm reinjecting cannula lure lock	1
0.7mm x 4cm reinjecting cannula lure lock	1
2.0mm x 5cm reinjecting cannula lure lock	1
1.2mm x 5cm reinjecting cannula lure lock	1
2.0mm x 15cm reinjecting cannula Super luer lock	1
3mm x 15cm V-Dissector Liposuction cannula lure lock	1
2.5mm x 20cm two lateral hole lure lock liposuction cannula	1
3.0mm x 20cm Blugerman I luer lock cannula	1
2.5mm X 20cm Harvester cannula with shevels hole lure lock	1
2.1mm x 20cm Tumescant Infiltrator luer lock cannula	1
3mm x 20cm Mercedes tip Super luer lock liposuction cannula	1
4mm x 20cm Spatula tip lure lock liposuction cannula	1
3mm x 25cm Tumescant Infiltrator luer lock cannula	1
4mm x 25cm Three standard hole lure lock liposuction cannula	1
Lure lock to Lure lock transfer adaptor	2
Lure lock to Lure lock transfer adaptor Steel	1
iv spike adapter	1
Handle for luer lock	1
Nano Fat with 3 filter	1
syringe Stopper	4
syringes Stop lock 10cc / 20cc	3

Nano Fat Transfer Set AD-5002

Specifically designed for nano fat and micro fat grafting procedures, this set provides everything needed to prepare and deliver refined fat with exceptional precision. The triple-filter nano fat processing unit ensures thorough emulsification of harvested tissue, while the colour-coded luer-to-luer transfer connectors in four sizes allow seamless, controlled transfer between syringes. A range of fine 1.2mm injector cannulas in varying lengths enables accurate placement of processed fat in superficial and delicate anatomical areas, making this set ideal for facial rejuvenation, skin quality improvement, and fine-detail contouring.

Description	Qty
Luer to luer , size 0.8mm 3 hole Geray	1
Luer to luer , size 1mm 3 hole light Blue	1
Luer to luer , size 1.5mm Green	1
Luer to luer , size 2mm Red	1
NANO FAT Transfer with 3 filter	1
Injector Cannula 1.2mm x 7cm	1
Injector Cannula 1.2mm x 10cm	1
Injector Cannula 1.0mm x 5cm	1
Injector Cannula 1.0mm x 7cm	1

All sets can be customized to suit your individual surgical needs. Cannulas and instruments may be added, removed, or replaced according to your preference. Please contact us for your personalized configuration.



Abdominal Cannula Set AD-5100



Purpose-built for abdominal liposuction, this set brings together a curated selection of high-performance cannulas to handle every stage of the procedure with confidence. Featuring a variety of specialized tip designs including Mercedes, Double Mercedes, Cobra K, Candy Cane, Spatula, and Becker Grater in diameters ranging from 3mm to 5mm and lengths up to 30cm, the set is engineered to tackle both superficial and deep abdominal fat with precision and efficiency. All cannulas are fitted with super handles for enhanced grip and control throughout the procedure.

Choose any Handle for this set



Description	Qty
Infiltrator Cannula 3mm x 15cm Super handle	1
Infiltrator Cannula 4mm x 30cm Super handle	1
Spatula Single port Cannula 3mm x 20cm Super handle	1
One hole standard Cannula 3mm x 20cm Super handle	1
Three diagonal hole Cannula 3mm x 30cm Super handle	1
Mercedes Liposuction Cannula 3mm x 30cm Super handle	1
Double Mercedes Cannula 5mm x 30cm Super handle	1
Cobra K Liposuction Cannula 4mm x 30cm Super handle	1
Candy Cane Cannula 4mm x 30cm Super handle	1
Becker grater round 5mm x 30cm Super handle	1

Full Body Liposuction Cannula Set AD-5300

A complete, all-encompassing set engineered to handle every stage of large-scale body contouring — from tumescent infiltration and fat harvesting to precision re-injection. The set features an extensive range of threaded and super luer cannulas in straight and curved configurations, including Mercedes, Basket, Toledo, Accelerator, and Tonnard designs, giving the surgeon the flexibility to adapt technique to tissue type and treatment zone. Fine single-hole fat injection cannulas and a nano fat filter support delicate re-injection and grafting, while luer-to-luer transfer adaptors, threaded connection handles, syringe stop locks, and an IV spike round out the kit for a fully self-contained, theatre-ready solution suited to comprehensive full body procedures.



Description	Qty	Description	Qty
Mercedes cannula 4x30cm Threaded Stright	1	Single hole fat injection Super luer 1.5x8cm	1
Mercedes cannula 4x30cm Threaded Curved	1	Single hole fat injection Super luer 2x7cm	1
Mercedes cannula 3x30cm Threaded Curved	1	Nano Fat Filter	1
Basket cannula 4x30cm Threaded Stright	1	Luer to Luer Transfer Adapter 2.4mm	1
Single Hole Fat injection 4x30cm Threaded Stright	1	Luer to Luer Transfer Adapter 2mm	1
Infiltration cannula 4x30cm Threaded Stright	1	Luer to Luer Transfer Adapter 1.2mm	1
Toledo cannula 4x30cm Threaded Stright	1	Luer to Luer Transfer Adapter 1mm	1
Accelerator cannula 4x30cm Threaded Stright	1	Sryinge Stop Lock 20cc	1
Mercedes Super luer 3x15cm curved	1	Luer Lock Sryinge Stoper Cap	4
Mercedes Super luer 2.5x15cm Curved	1	Female luer needle conector	1
Infiltration Super luer 2.5x15cm	1	Adaptor handle with connection threaded	1
Infiltration Super luer 2.5x15cm Curved	1	Adaptor handle with hole, for interruption of suction	1
Tonnard Super luer 2.5x20cm	2	Adaptor with luer lock connection	1
Single hole fat injection Super luer 1.2x8cm	1	Iv spike	1
Large Connection handle for Machine suction	1		

Breast Contouring Threaded Cannula Set AD-5303



Description	Qty
Single hole 3x20cm threaded Curved	1
Single hole 3x20cm threaded Curved	1
Single hole 3x25cm threaded Curved	1
Single hole 3x25cm threaded Stright	1
Three hole cannula 3x20cm threaded	1
V Dissector Cannula 3x20cm threaded	1
Infiltration Cannula 3x25cm threaded	1
One Central two lateral hole 3x25cm threaded	1
Adaptor handle with connection threaded	1
Female luer needle conector	1

A dedicated set designed for precision liposuction and fat grafting of the breast, built around an internally threaded cannula design that offers a secure and stable connection between cannula and handle for enhanced control. The set includes a range of 3mm threaded cannulas in straight and curved configurations, including single-hole, three-hole, V-Dissector, infiltration, and combination central-lateral hole designs — paired with a matching threaded connection handle for quick and reliable attachment. A female luer needle connector is also included, rounding out the set for refined breast contouring procedures.

Abdominal Threaded Cannula Set AD-5304

Description	Qty
Mercedes cannula 4x30cm Threaded Stright	1
Mercedes cannula 4x30cm Threaded Curved	1
Mercedes cannula 3x30cm Threaded Curved	1
Basket cannula 4x30cm Threaded Stright	1
Single Hole Fat injection 4x30cm Threaded Stright	1
Infiltration cannula 4x30cm Threaded Stright	1
Toledo cannula 4x30cm Threaded Stright	1
Accelerator cannula 4x30cm Threaded Stright	1
Female luer needle conector	1
Adaptor handle with connection threaded	1
Adaptor handle with hole, for interruption of suction	1

Engineered for high-volume abdominal liposuction, this set is built around a secure internally threaded cannula design that ensures a stable, reliable connection between cannula and handle throughout the procedure. The set features a strong selection of 4mm and 3mm threaded cannulas in straight and curved configurations, including Mercedes, Basket, Toledo, Accelerator, Infiltration, and Single Hole Fat Injection designs, to handle both tunneling and aspiration of dense abdominal fat with efficiency. Two specialized adaptor handles, including one with a suction-interruption hole, along with a female luer needle connector, complete the set for a versatile and controlled abdominal contouring experience.





Address:

Adas Instruments

Nawan Pind, P.O Chitti Sheikhan

Sialkot-51310,

Pakistan.

Mob. : +92 321 6126225

Email: info@adasinstrumentsonline.com adas.surgical@gmail.com

www.adasinstrumentsonline.com